

Atomic Physics DAQ Container


- Phone: 2938
- PC: SDPC027 (140.181.77.246)
- PC: SDPD065, Optiplex 755 (140.181.115.078)
MAC: 00219B4D7D5E
- Network sockets: TH2B 01D19 (SDPC027)
TH2B 01D20 (SDPC065)
- Power sockets: 2x, label: = PE3
(left wall next to door) >IL04.6-1F3
- Contact person:
Angela Bräuning-Demian, 2711

Racks 1 + 2




- Rack 2 supplied by rack 3 !!!!!!!
- Be careful when switching off rack 3!!!!!!

Racks 3 + 4 (TOP missing)




- Two 2-pin Lemo cables for timing from SIS timing module in other electronics room of container

Racks 3 (TOP)




- Cables connect to ????

Racks 4 (Patch panel)



- 8x RG216
- 4x SHV
- 2x 25-pol.

Racks 4 (Network I)



- Frau Vincelli in charge???
- Connection to ????
- Optical link to SD electronics room in BG2 ??

Racks 4 (Network II)

